

# IP-2030

## Vital sign monitor

Real-time measurement patient vital signs



Advanced SpO2 technology  
(Motion artifact)



Optional NIBP / temperature  
/ ECG / EtCO2



5 inch high resolution  
screen



Compact & light weight




Mounting option  
(IV pole clamp)



Perfusion range:  
0.05%-20%

Intuitive approach to patient vital signs measurement.  
Basic pulse oximeter & patient monitoring use



Type A : SpO2 + TEMP




Type B : SpO2 + TEMP



Type C : SpO2 + NIBP + TEMP + ECG



Type D : SpO2 + NIBP + ECG + Respiration + EtCO2



**ECG**

- Lead: 3 or 5 leads
- Heart rate range 20bpm to 300bpm

**SpO2**

- Pulse rate: 25 to 250bpm

**NIBP**

- Maximized sensitivity of NIBP
- The accurate measurement from adult to neonate

## ECG Analysis

PVC and detection of 12 kinds of arrhythmia

## SpO2

Perfusion index: 0.05~20%

## EtCO2 Module

Built in sidestream / mainstream EtCO2 module

Interface: USB, nurse call, wi-fi (option)

Suitable for ambulatory use (fast booting time)



## Vital sign monitor (Can be used alone)

Suitable for spot check & emergency use

Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800  
F +82-2-546-8346  
E sales@infunix.com

Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500  
E genoray@genoray.com  
W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729  
E inquiry@genorayamerica.com  
W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498  
E info@genorayeuropa.com  
W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971  
E Info@genorayjapan.co.jp  
W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37  
+90-216-463-57-53  
E genorayturkiye@genoray.com

GENORAY Shanghai

T +86-021-6785-8538  
W www.genoray.com.cn