



IP-4050

PureScope series

Bedside monitor



7 Parameter
6 waveforms



10.4 inch
touch screen



ECG (3 or 5 lead), SpO2, NIBP,
2-IBP, 2-temperature, respiration
(Opt.EtCO2 / printer)





Detects 12 types of
ECG arrhythmia and
displays PVC and ST level



Notification system when
it detects an abnormality
in a patient



Central monitoring system
(Monitoring 32 beds)



ECG analysis

- 12 Kinds of arrhythmia
- Arrhythmia alarm
- Display of PVC & ST level
- Pacemaker detection

Trend

72 hours of graphic and tabular trend

Built-in EtCO2 module

- Easy connect EtCO2 cannula to the monitor directly
- Convenient & user-friendly use
- Prevent from destroy or loss of external module

Optional

- Touch screen
- Rolling stand & wall mount
- Printer



Central monitoring system

Bidirectional communication with patient monitors
Monitoring 4, 8, 16, 32 beds
Waveform, parameter and alarm are saved.
Support hl7

| Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800
F +82-2-546-8346
E sales@infunix.com

| Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500
E genoray@genoray.com
W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729
E inquiry@genorayamerica.com
W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498
E info@genorayeuropa.com
W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971
E Info@genorayjapan.co.jp
W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37
+90-216-463-57-53
E genorayturkiye@genoray.com

GENORAY Shanghai

T +86-021-6785-8538
W www.genoray.com.cn