


IP-5000

PureScope series

Multi parameter patient monitor



7 Parameter
9 waveforms



15 Inch
touch screen



Configurable color of
parameters & waveforms





Detects 12 types of
ECG arrhythmia and
displays PVC and st level



Operates 4 hours
charged battery



Central monitoring system
(Monitoring 32 beds)



ECG analysis

- 12 Kinds of arrhythmia
- Arrhythmia alarm
- Display of PVC & ST level
- Pacemaker detection

Trend

168 Hours of graphic and tabular trend

SpO2 averaging (4/8/12/16)

Alarm duration adjustable

Optional

- Touch screen
- Rolling stand & wall mount
- Printer



Central monitoring system

Bidirectional communication with patient monitors
Monitoring 4, 8, 16, 32 beds
Waveform, parameter and alarm are saved.
Support hl7

| Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800
F +82-2-546-8346
E sales@infunix.com

| Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500
E genoray@genoray.com
W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729
E inquiry@genorayamerica.com
W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498
E info@genorayeuropa.com
W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971
E Info@genorayjapan.co.jp
W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37
+90-216-463-57-53
E genorayturkiye@genoray.com

GENORAY Shanghai

T +86-021-6785-8538
W www.genoray.com.cn